

# OptiVac® AC/DC PORTABLE SUCTION UNIT

## Model G180



Model G180 with Disposable Hydrophobic Bacteria Filter and optional Carrying Case

**THE QUALITY OF GOMCO® WHERE YOU NEED IT MOST** – The G180 is the second generation OptiVac in the industry-leading Gomco line of portable suction equipment. Its compact design, ease of operation, rechargeable battery, AC power supply, and fully-enclosed unit are just a few of its advanced features.

### FEATURES OF THE G180 MODEL

- Run time operation on fully-charged battery for approximately 3 hours at maximum vacuum.
- 25-550 mm Hg vacuum range, coupled with greater than 30 lpm flow, provides powerful suction for faster aspiration.
- Fully chargeable in 8 hours.
- Will run on AC power even with low battery.
- Will run for approximately 1 hour after charged for just 2 hours.
- Long-life rechargeable battery eliminates the need for frequent replacement.
- Sealed water-resistant unit is easy to clean, and minimizes potential contamination or damage to the compressor and battery.
- Control panel at an angle for easy viewing.
- Quality-engineered thermoplastic provides impact resistance and durability.
- 3-year warranty on the unit, and 1 year on the battery.

MODEL	COLLECTION BOTTLE	MOTOR	COMPRESSOR	VACUUM RANGE	FLOW RATE	DIMENSIONS HxDxW	SHIPPING WEIGHT
G180	1500 ml Disposable Collection Canister	3.4 amps	Diaphragm	1-22 in. Hg 25-550 mm Hg	30 lpm minimum	9.4 x 7.5 x 16.8 in. 23.9 x 19.1 x 42.7 cm	11.4 lbs. 5 kg
G180-Jar	Half-gallon Glass Bottle	3.4 amps	Diaphragm	1-22 in. Hg 25-550 mm Hg	30 lpm minimum	9.4 x 7.5 x 16.8 in. 23.9 x 19.1 x 42.7 cm	11.4 lbs. 5 kg

### ORDERING INFORMATION

G180	OptiVac Portable Suction Unit with Gauge and Regulator, ON & OFF Buttons, one plastic 1500 ml disposable collection canister with lid, one length of 15" clear flexible tubing (pump tube), one length of 72" clear flexible tubing (patient tube), three hydrophobic bacteria filters, AC power cord with hospital-grade plug; 115/60 voltage
G180CE	Same as G180; 230/50 voltage
SU1060	Battery for G180 (12 volt DC, 5 amp-hour)
L200023	AC Power Cord
L190-BAG	Carrying Case with shoulder strap
20-08-0003	Disposable Collection Canister, 1500 ml, stem inlet, 48/ case
20-08-0016	Disposable Collection Canister, 1500 ml, stem inlet, 16/case
01-90-3928	Disposable Hydrophobic Bacteria Filter, 1/4"x1/4" to barb, 3/pkg.
S610100	Suction Tubing Kit (one each 72" and 15" length, 1/4" ID, clear PVC)

### ORDERING INFORMATION - MODEL 305

w/Glass Bottle	w/Disposable Canister	Voltage
01-22-0305*	01-22-3005*	115/60

\*Tested to 60601-1 IEC Standard; AC power only model - no battery