

# THORACIC DRAINAGE PUMPS

Models 6020, 6050



**THORACIC DRAINAGE PUMPS** – Gomco® mobile units are designed as a special aid to the chest surgeon for post-operative care.

Stands include 3" non-marking casters, handy cord storage, chemical resistant stainless steel top, compact stand design, graduated 2800 ml polycarbonate plastic collection bottle, and hospital grade plug. Models include illuminated on/off switches.

**MODEL 6020** – This model is designed for post-operative care of chest wounds. The manometer is designed to give a visual indication of function, and to provide regulation of suction from 0 to 25 cm of water while operating for extended periods of time. This pump has applications for neonatal care

or for situations where low volume and flow rates are required. Single 2800 ml collection bottle with calibrated water seal tube is included. This unit is UL listed.

**MODEL 6050** – This unit provides regulation of suction (drainage) from 0 to 50 cm of water for rapid removal of fluids and gas from the pleural cavity. This pump offers high flow rates and a dual 2800 ml polycarbonate plastic collection bottle system with two calibrated water seal tubes and two drainage tubes. A cap is provided when only one drainage tube is required. This unit is UL listed.

MODEL	COLLECTION BOTTLE	COMPRESSOR	VACUUM RANGE	FLOW RATE	DIMENSIONS HxDxW	SHIPPING WEIGHT
6020	Polycarbonate Plastic (one)	Thermotic Pump	0-25 cm H <sub>2</sub> O	2.3 lpm (minimum)	33 x 17 x 18 in. 84 x 43 x 46 cm	65 lbs. 29 kg
6050	Polycarbonate Plastic (two)	Continuous Duty High-Volume Pump	0-50 cm H <sub>2</sub> O	110 lpm (minimum)	33 x 17 x 18 in. 84 x 43 x 46 cm	63 lbs. 28 kg

## ORDERING INFORMATION

w/Polycarbonate

Plastic Bottle	Voltage
01-31-6020	110/50
01-32-6020	115/60
01-34-6020	220/60
01-73-6020	220/50
01-31-6050	110/50
01-32-6050	115/60
01-33-6050	220/50
01-34-6050	220/60